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HTM 01-01

Position of PQ sensors

2.40 Temperature sensors should be as described in HTM 01-01 Part B.

2.41 Temperature sensors should be placed in the following positions:

- a. one on/in each of three items that are slowest to attain the sterilization temperature;
- b. one on/in each of three items that are fastest to attain the sterilization temperature;
- c. if the load consists of fewer than six items, one on/in each item;
- d. if the load includes lumen devices, temperature sensors should be placed to monitor the environment within the lumen rather than the device's surface, at the

most challenging position within the device. In cases where temperature cannot be used to determine the presence of residual air (for example a narrow lumen or metal device in which the residual air rapidly attains steam temperature), alternative sensor technology should be used. Examples include chemical and biological indicators.

HTM 01-01

Microbiological test for PQ

2.51 This test is designed to be used in exceptional circumstances as an additional PQ test for steam sterilizers. The microbiological test should ideally follow a satisfactory thermometric test, using the identical loading condition and operating cycle. There may be situations where thermometric tests are not

possible, for example with narrow-lumened instruments, where it is not physically possible to place a thermocouple or temperature sensor into the lumen without altering the nature of the load. Reference should be made to BS EN 556-1 for sterility assurance requirements.

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Technical and project manager
09-68241222
michael.paul@getinge.com

www.getinge.com

Getinge is a leading global provider of innovative solutions for operating rooms, intensive-care units, hospital wards, sterilization departments and for life science companies and institutions. Based on our first-hand experience and close partnerships, we are improving the every-day life for people, today and tomorrow.

GETINGE

PASSION FOR LIFE